


The Steripharm Face Mask is intended for use by members of the general public as well as HCPs in healthcare settings as a face mask to cover the nose and mouth for source control (to help contain the wearer's respiratory secretions). It may be worn and safely discarded after a single use.



- Optimizing Safety
- Lightweight
- Comfortable for long durations
- Soft ear loops for easier donning
  - Non-irritating materials used for manufacturing
  - Body Contacting Materials
  - Nonwoven Polypropylene
  - Cotton
  - Plastic Coated Steelwire
- Compliant with Global Standards



## Warnings

- Do Not Use in clinical settings with significant risk of infection exposure through inhalation.
- Do Not Use in surgical environments
- Do Not Use in environments where risk to hazardous fluid exposure is greater.
- Do Not Use in the presence of highly flammable sources.
- Do Not Reuse.

### Steripharm Ltd. FZC

Plot No. L314, Fujairah Free Zone

Fujairah, United Arab Emirates

Telephone: Tel : 00971 9 2228862

Website: [www.steripharmuae.com](http://www.steripharmuae.com)

Email: [customersupport@steripharm.ae](mailto:customersupport@steripharm.ae)

The product has not been FDA cleared or approved. The product has been authorized by FDA under an EUA for use as source control by the general public as well as by HCP in healthcare settings. The product is authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of medical devices, including alternative products used as medical devices, during the COVID-19 outbreak, under section 564(b)(1) of the act, 21 U.S.C. 360bbb3(b)(1) unless the authorization is terminated or revoked sooner."